REF 5800102	Technofluor ADAMTS13 Activity 0 CAL 2 x 1 mL	ENG
		DEU
		FRA
		NED
		ITA
		ESP
3050526A Rev.003 10/10/2023		

Technofluor ADAMTS13 Activity 0 CAL - English
INTENDED USE
Calibration of the Technofluor ADAMTS13 activity assay on Ceveron s100 in combination with Coagulation Reference.
SUMMARY

Von Willebrand factor (VWF) promotes clot formation by tethering platelets at the site of vessel injury and can also contribute to platelet aggregation. It is comprised of subunits of differing sizes, termed multimers, whose adhesiveness increases with size. The enzyme ADAMTS13 (a disintegrin-like and metalloproteinase with thrombospondin type 1 motif 13) acts as a gatekeeper against generation of highly thrombogenic ultra large multimers by regulating multimer size via a specific cleavage site, thereby breaking down large multimers into smaller, less reactive forms. Ultra large multimers can achieve millimeters in length if unregulated. In ADAMTS13 deficiency states, unusually large hyperfunctional VWF multimers can accumulate, leading to microvascular thrombosis and organ damage/failure due to VWF-platelet aggregates. In turn, a microangiopathic haemolytic anaemia can ensue, red blood cells being sheared as they travel past and through the microthrombi. ADAMTS13 deficiency can be hereditary or acquired, the latter commonly due to ADAMTS13 autoantibodies, and is termed thrombotic thrombocytopenic purpura.

In the Technofluor assay, ADAMTS13 from the plasma sample cleaves a VWF73 based substrate, thereby releasing the dark quencher. Consequently, the energy transfer that quenches fluorescence from an adjacent emitter does not occur, allowing emission of fluorescence. The fluorescent signal is proportional to the ADAMTS13 activity in the sample.

REAGENTS		
The Technofluor ADAMTS13 Activity 0 CAL contains:		
Reagent / Content	Description	
2 x 1 mL	0 calibrator (CAL)	Lyophilized human plasma with no ADAMTS13 activity

Material required (not supplied with the kit)

- Distilled water
- Precision pipettes
- Test systems
- REF 5800100 Technofluor ADAMTS13 Activity

 25 tests
- Calibration Plasma*
- REF 5220110 Coagulation Reference

 5 x 1 mL
- Laboratory timer

* or any other package sizes.

Warning and precautions

- IVD for in vitro diagnostic use only.
- This kit is intended for use by personnel trained in laboratory procedures and universal precautions for the use of chemicals and potentially biohazardous substances must be applied.
- All human blood or plasma products as well as test samples must be considered as potentially infectious. They have to be handled with appropriate care and in strict observance of safety regulations. The rules pertaining to disposal are the same as applied to disposing hospital waste.
- Calibrators and control plasmas are made from human blood and any individual plasma involved in the procedure is tested HbsAg, HIV 1/2 Ab and HCV-Ab-negative. However, all human blood products should be handled as potentially infectious material.
- Get a Material Safety Data Sheet for this product from www.technoclone.com.

Stability and storage

The expiry date printed on the labels is only applicable to storage of the unopened containers at 2...8 °C.

Reagent / Content	Ceveron s100 (open vial)	2...8 °C (closed vial)	< -20 °C (closed vial)
0 calibrator (CAL)	8 hours	-	2 months

Calibrator should only be frozen once. Thawing must be performed rapidly in a water bath kept at 37 °C.

TEST PROCEDURE
Preparation of reagents

Before starting the test, all the required components are to be brought to room temperature.

Avoid foam formation when reconstituting plasmas and mixing reagents or buffers.

- Calibrator*: Dissolve each bottle of lyophilized calibrator in 1.0 mL distilled water and swirl gently. Allow the reconstituted material to stand for 10 minutes at room temperature before use. For standardizing tests, a reconstitution time of 30 minutes is recommended. Swirl to mix before use.

Performance of the test
The Technofluor ADAMTS13 Activity is performed on Ceveron s100 with the respective application.
Calibration is performed using a serial dilution of Coagulation Reference in Technofluor ADAMTS13 Activity 0 CAL. Normal and abnormal controls are recommended for a complete quality control program. Coagulation Control N and Coagulation Control A are designed for this program. Each laboratory should establish its own mean and standard deviation for a quality control program in order to monitor laboratory testing. Controls should be analyzed before validating patient results in accordance with good laboratory practice.

STANDARDISATION
Technofluor ADAMTS13 Activity 0 CAL and Coagulation Reference (Calibrator) are calibrated using Technofluor ADAMTS13 Activity and traceable to the WHO International Standard for ADAMTS13 Activity. Consult the respective batch table for each material.
LITERATURE
Please contact Technoclone www.technoclone.com or your local distributor.
EDITORIAL NOTE

This document is available in several languages. The translations have been done using the master document in English. In the event of doubts or discrepancies, the wording in the master document in English shall take precedence.
Technofluor ADAMTS13 Activity 0 CAL - Deutsch
ANWENDUNG
Das Technofluor ADAMTS13 Activity 0 CAL dient in Kombination mit dem Coagulation Reference der Kalibration des Technofluor ADAMTS13 Aktivitätskit am Ceveron s100.
ZUSAMMENFASSUNG

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Technofluor ADAMTS13 Activity 0 CAL - Français

Technofluor ADAMTS13 Activity 0 CAL - German

Technofluor ADAMTS13 Activity 0 CAL - Spanish

Technofluor ADAMTS13 Activity 0 CAL - Português

Technofluor ADAMTS13 Activity 0 CAL - Italiano

Technofluor ADAMTS13 Activity 0 CAL - Nederlands

Technofluor ADAMTS13 Activity 0 CAL - Norsk

Technofluor ADAMTS13 Activity 0 CAL - Svenska

Technofluor ADAMTS13 Activity 0 CAL - Slovenian

Technofluor ADAMTS13 Activity 0 CAL - Esperanto

Technofluor ADAMTS13 Activity 0 CAL - Chinese

Technofluor ADAMTS13 Activity 0 CAL - Japanese

Technofluor ADAMTS13 Activity 0 CAL - Korean

Technofluor ADAMTS13 Activity 0 CAL - Russian

Technofluor ADAMTS13 Activity 0 CAL - Polish

Technofluor ADAMTS13 Activity 0 CAL - Czech

Technofluor ADAMTS13 Activity 0 CAL - Hungarian

Technofluor ADAMTS13 Activity 0 CAL - Greek

Technofluor ADAMTS13 Activity 0 CAL - Hebrew

Technofluor ADAMTS13 Activity 0 CAL - Hindi

Technofluor ADAMTS13 Activity 0 CAL - Malayalam

Technofluor ADAMTS13 Activity 0 CAL - Marathi

Technofluor ADAMTS13 Activity 0 CAL - Nepali

Technofluor ADAMTS13 Activity 0 CAL - Punjabi

Technofluor ADAMTS13 Activity 0 CAL - Sinhala

Technofluor ADAMTS13 Activity 0 CAL - Tamil

Technofluor ADAMTS13 Activity 0 CAL - Telugu

Technofluor ADAMTS13 Activity 0 CAL - Urdu

Technofluor ADAMTS13 Activity 0 CAL - Vietnamese

Technofluor ADAMTS13 Activity 0 CAL - Bengali

Technofluor ADAMTS13 Activity 0 CAL - Gujarati

Technofluor ADAMTS13 Activity 0 CAL - Kannada

Technofluor ADAMTS13 Activity 0 CAL - Malay

Technofluor ADAMTS13 Activity 0 CAL - Odia

Technofluor ADAMTS13 Activity 0 CAL - Oriya

Technofluor ADAMTS13 Activity 0 CAL - Sinhala

Technofluor ADAMTS13 Activity 0 CAL - Thai

Technofluor ADAMTS13 Activity 0 CAL - Vietnamese

Technofluor ADAMTS13 Activity 0 CAL - Chinese

Technofluor ADAMTS13 Activity 0 CAL - Japanese

Technofluor ADAMTS13 Activity 0 CAL - Korean

Technofluor ADAMTS13 Activity 0 CAL - Russian

Technofluor ADAMTS13 Activity 0 CAL - Polish

Technofluor ADAMTS13 Activity 0 CAL - Czech

Technofluor ADAMTS13 Activity 0 CAL - Hungarian

Technofluor ADAMTS13 Activity 0 CAL - Greek

Technofluor ADAMTS13 Activity 0 CAL - Hebrew

Technofluor ADAMTS13 Activity 0 CAL - Hindi

Technofluor ADAMTS13 Activity 0 CAL - Malayalam

Technofluor ADAMTS13 Activity 0 CAL - Marathi

Technofluor ADAMTS13 Activity 0 CAL - Nepali

Technofluor ADAMTS13 Activity 0 CAL - Punjabi

Technofluor ADAMTS13 Activity 0 CAL - Sinhala

Technofluor ADAMTS13 Activity 0 CAL - Tamil

Technofluor ADAMTS13 Activity 0 CAL - Telugu

Technofluor ADAMTS13 Activity 0 CAL - Urdu

Technofluor ADAMTS13 Activity 0 CAL - Vietnamese

Technofluor ADAMTS13 Activity 0 CAL - Bengali

Technofluor ADAMTS13 Activity 0 CAL - Gujarati

Technofluor ADAMTS13 Activity 0 CAL - Kannada

Technofluor ADAMTS13 Activity 0 CAL - Malay

Technofluor ADAMTS13 Activity 0 CAL - Odia

Technofluor ADAMTS13 Activity 0 CAL - Oriya

Technofluor ADAMTS13 Activity 0 CAL - Sinhala

Technofluor ADAMTS13 Activity 0 CAL - Thai

Technofluor ADAMTS13 Activity 0 CAL - Vietnamese

Technofluor ADAMTS13 Activity 0 CAL - Chinese

Technofluor ADAMTS13 Activity 0 CAL - Japanese

Technofluor ADAMTS13 Activity 0 CAL - Korean

Technofluor ADAMTS13 Activity 0 CAL - Russian

Technofluor ADAMTS13 Activity 0 CAL - Polish

Technofluor ADAMTS13 Activity 0 CAL - Czech

Technofluor ADAMTS13 Activity 0 CAL - Hungarian

Technofluor ADAMTS13 Activity 0 CAL - Greek

Technofluor ADAMTS13 Activity 0 CAL - Hebrew

Technofluor ADAMTS13 Activity 0 CAL - Hindi

Technofluor ADAMTS13 Activity 0 CAL - Malayalam

Technofluor ADAMTS13 Activity 0 CAL - Marathi

Technofluor ADAMTS13 Activity 0 CAL - Nepali

Technofluor ADAMTS13 Activity 0 CAL - Punjabi

Technofluor ADAMTS13 Activity 0 CAL - Sinhala

Technofluor ADAMTS13 Activity 0 CAL - Tamil

Technofluor ADAMTS13 Activity 0 CAL - Telugu

Technofluor ADAMTS13 Activity 0 CAL - Urdu

Technofluor ADAMTS13 Activity 0 CAL - Vietnamese

Technofluor ADAMTS13 Activity 0 CAL - Bengali

Technofluor ADAMTS13 Activity 0 CAL - Gujarati

Technofluor ADAMTS13 Activity 0 CAL - Kannada

Technofluor ADAMTS13 Activity 0 CAL - Malay

Technofluor ADAMTS13 Activity 0 CAL - Odia

Technofluor ADAMTS13 Activity 0 CAL - Oriya

Technofluor ADAMTS13 Activity 0 CAL - Sinhala

Technofluor ADAMTS13 Activity 0 CAL - Thai

Technofluor ADAMTS13 Activity 0 CAL - Vietnamese

Technofluor ADAMTS13 Activity 0 CAL - Chinese

Technofluor ADAMTS13 Activity 0 CAL - Japanese

Technofluor ADAMTS13 Activity 0 CAL - Korean

Technofluor ADAMTS13 Activity 0 CAL - Russian

Technofluor ADAMTS13 Activity 0 CAL - Polish

Technofluor ADAMTS13 Activity 0 CAL - Czech

Technofluor ADAMTS13 Activity 0 CAL - Hungarian

Technofluor ADAMTS13 Activity 0 CAL - Greek

Technofluor ADAMTS13 Activity 0 CAL - Hebrew

Technofluor ADAMTS13 Activity 0 CAL - Hindi

Technofluor ADAMTS13 Activity 0 CAL - Malayalam

Technofluor ADAMTS13 Activity 0 CAL - Marathi

Technofluor ADAMTS13 Activity 0 CAL - Nepali

Technofluor ADAMTS1

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Ceveron is a registered trademark of Technoclone