# Siron LS (aPTT liquid)

GB

IT

PT

DE

ES

FR

REF	5035105	Siron LS (aPTT liquid)	2 x 4 mL
REF	5035107	Siron LS (aPTT liquid)	10 x 4 mL
REF	5035109	Siron LS (aPTT liquid)	10 x 10 mL

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Ključova slova / Značenje simbola			
	manufacturer / Hersteller / fabbricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produzent / Κατασκευαστής / Производитель / Производител / výrobce / Proizvođač		expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace/ срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / oppbevaringstemperatur / Oppbevaringstemperatur / θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg brugsvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potřeba řídit se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj upustvo pre upotrebe
	CE-mark / CE-Kennzeichnung / marchio CE / marca de CE / Simbolo da CE / marquage CE / CE-mærkning / CE-mærket / CE-merke / CE-σημάδι / CE марка / CE-označení / маркировка CE / značka CE / CE-marca		determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelse / Bestemmelse / προσδιορισμοί / брой тестове / stanovení / определний / роčet stanovení / Definicija
<b>AQUA</b>	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / децилирана вода / destilovaná voda / дистиллированная вода / destilovaná voda / Destilana Voda	<b>LOT</b>	lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / партида номер / šarže / лот / šarže / Serija
<b>BUF</b>	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionspuffer / Reaktionspuffer / Reaktionspuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakční pufr / Reakcioni pufer	<b>MTP</b>	microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροτιτλοδότησης / Микротитрна плака / Микропланшет / Mikrotitrační destička / Mikrotitracione ploče
<b>CAL</b>	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	<b>REF</b>	catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
<b>CONJ</b>	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat	<b>RTU</b>	ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk / έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
<b>CONT</b>	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola	<b>STOP</b>	stop solution / Stoppløsning / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stoppløsning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп раствор / Стоп-раствор / Zastavovací roztok / Stop solucija
<b>DIL</b>	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / spädd eller upplös i / fortyndes eller opløses i / Fortyndes eller oppløses i / αραιωση ή διάλυση σε / растворяете или разредете с / zředit anebo rozpustit v / разбавить или растворить в / nafedte nebo rozpustite v / razrediti ili rastvoriti u	<b>SUB</b>	substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
<b>INC</b>	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationspuffer/ Inkubationspuffer/ Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer	<b>WASH</b>	washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningsskoncentrat / Vaskeopløsningskoncentrat / vaskeløsningskoncentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промывочного раствора / Koncentrát promývacieho roztoku / Koncentrat solucije za ispiranje
<b>IVD</b>	in vitro diagnostic use / in vitro Diagnostikum / diagnostico in vitro / diagnóstico en vitro / diagnóstico in vitro / diagnostic in vitro / för in vitro diagnostik / in vitro diagnostik / in vitro diagnostisk bruk / рηση διαγνωστικής εντόβ σωλήνα / за ин vitro диагностика / pro in vitro diagnostiku / использовать для диагностики in vitro / diagnostický prostředek in vitro / in vitro dijagnostika		

**PRODUCT DESCRIPTION****INTENDED USE**

Siron LS is an activated partial thromboplastin time (aPTT) reagent to be used:

- as a screening test for the intrinsic coagulation system
- in the diagnosis and treatment of haemophilia A and B
- for specific determinations of factors VIII, IX, X, XI and XII when used in combination with the corresponding deficiency plasma, for control of heparin treatment
- for detection of Lupus anticoagulant

Siron LS measures the overall activity of the intrinsic coagulation system. Siron LS is sensitive to the factors XII, XI, X, IX, VIII, V and II, but insensitive to factors VII and XIII.

**COMPOSITION**

Siron LS is an activated partial thromboplastin time (aPTT) test. This liquid preparation is a stable aqueous suspension of phospholipids. FXII activation of Siron LS is achieved by addition of ellagic acid.

**MATERIAL REQUIRED** (not supplied with the kit)

- Buffer\*\*

[REF] 5277015 CaCl<sub>2</sub> 25 mmol/L 100 mL

- Control Plasma Normal and Abnormal \*\*

[REF] 5020040 Coagulation Control N 5 x 1 mL

[REF] 5021055 Coagulation Control A 5 x 1 mL

[REF] 5220130 Coagulation Reference 5 x 1 mL

\*\* or any other package sizes, special Ceveron® alpha or TECHNOCLOT® Control and Calibration reagents of Technoclone.

**WARNING AND PRECAUTIONS**

- IVD for *in vitro* diagnostic use
- Siron LS contains no human material. FXII activation of Siron LS is achieved by ellagic acid and phospholipids from synthetic origins.
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.

**STABILITY AND STORAGE**

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8°C.

Stability opened vials:

+2...8°C	+12°C (Ceveron® alpha)	RT*
4 weeks	10 days	24 hours

Upon storage, caps should be screwed tightly.

\*=room temperature (+18...25°C)

**TEST PROCEDURE****PREPARATION OF PLASMA SAMPLES**

Plasma separation: Mix 9 parts of venous blood and 1 part of Sodium Citrate Solution (0.11 mol/L or 0.13 mol/L) and centrifuge for 15 minutes at a RCF of at least 2500 g (corresponding to DIN 58905) to obtain platelet poor plasma.

**STORAGE OF PLASMA**

Samples of plasma are stored in plastic tubes or siliconized vials at room temperature before testing:

- For patients that are not on heparin therapy, within **3 hours**
- For patients on heparin therapy, within **1 hour**
- For longer period - 20°C, 2 weeks  
- 70°C, 6 months

**PREPARATION OF REAGENT**

Each vial contains 4 mL / 10 mL ready to use aPTT reagent. Before each use, let the reagent stabilize for 30 min at room temperature (+18...25°C) while shaking the vial from time to time.

**PERFORMANCE OF THE TEST****CEVERON**

Technoclone provides application sheets for Ceveron® alpha. The application sheets contain analyser/assay specific handling and performance information which may differ from that provided in this instruction for use. In this case the information contained in the application sheets supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

**MANUAL**

The CaCl<sub>2</sub> 25 mmol/L should be preheated to +37°C.

Pipetting Scheme:

100 µL plasma sample
+ 100 µL Siron LS in a pre-warmed tube
incubate for <b>3 minutes</b>
+ 100 µL CaCl <sub>2</sub> solution (25 mmol/L; +37°C) Start stop-watch simultaneously
Leave the tube for 20 sec. before tilting. Determine the clotting end point when clot is formed. The results are given in seconds.

**ANALYSES RESULTS****REFERENCE RANGE**

**Normal range:** 25 - 38 seconds

The normal range refers to analyzers with optical detection (Ceveron® alpha). It is recommended that each laboratory establish and controls its normal range.

**Siron LS in control of heparin therapy:**

Siron LS seconds of 2 - 4 times the normal has been found to correspond to the accepted therapeutic level of heparin. However, the actual therapeutic level should be decided from the clinical situation.

**STANDARDISATION**

No international calibrator is available for the standardization of aPTT tests. Siron LS documentation is based on studies testing plasma samples from normal subjects and heparinised plasma samples. An in-house reference batch is established in order to avoid batch-to batch variation.

**LIMITATION OF THE TEST**

- The blood samples should be taken in non-activating samples tubes in order to avoid activation of the clotting factors in the plasma sample.
- Plasma samples should be stored at room temperature to avoid cold activation.
- Platelet poor plasma samples should be used.
- It is important to use the correct incubation time, which is 3 minutes for Siron LS when testing both heparinised and normal plasma samples.
- The calcium concentration used for re-calcification should be 25 mmol/L when testing normal, abnormal or heparin plasmas.
- The REFERENCE RANGE represents approximate values only. Each laboratory should determine its own normal range to take account of laboratory specific factors, such as the type of coagulometer used, reagent temperature, patient population, etc...
- The values found when testing control plasma normal and control plasma abnormal should be compared to the value given on the data-key for the corresponding lot of Control plasma.
- If the results obtained are outside the recommended range, avoid measuring patient samples until the problem is solved.
- A new calibration is required for each batch of Siron LS and additional reagents (e.g. deficient plasmas) where a calibration curve is necessary and for each instrument used. Also a new calibration is recommended, if software changes are introduced or following a major service of either instruments or equipment.
- Clotting time is not affected in this assay by samples containing ≤ 1000 mg/dL Haemoglobin, ≤46 mg/dL unconjugated bilirubin and ≤ 1000 mg/dL triglycerides. However elevated levels of conjugated bilirubin causes a prolonged clotting time.

**PERFORMANCE CHARACTERISTICS**

Performance data are given below. Results obtained in individual laboratories may differ.

**PRECISION**

Reproducibility was determined with different samples (in series and day to day). The following results were obtained:

Sample	Intra assay		Inter assay	
	Sample 1	Sample 2	Sample 1	Sample 2
n	21	21	20	20
Mean (sec)	28.2	73.3	28.2	65.7
CV (%)	0.99	0.92	1.96	2.3

**COMPARISON OF METHODS OR CORRELATION**

Following correlation was obtained in comparing SIRON LS:

$$\text{Actin® FSL} \quad y = 0.775x + 7.9 \quad R^2 = 0.95$$

**LITERATURE**

Please contact Technoclone or your local distributor.

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