REF 5010024 AK-Verification Kit

3015102 Rev. 001 [01] 07/2015

AK-Verification Kit - English

INTENDED USE

The AK-Verification Kit is a tri-level quality control set for the verification of International Normalized Ratio (INR) reporting.

SUMMARY

Verification of INR results in accordance to international guidelines^{1,2,3,4} requires a minimum of three certified plasmas within an INR range from 1.5 to 4.5. Each time a new ISI is calculated or a new PT/INR calibration curve is generated, the system should be verified before patient INR values are determined.

REAGENTS

The AK-Verification Kit contains three plasma levels (1,2, and 3) which are produced exclusively from the plasma of donors on long-term oral anticoagulant treatment. As a result, the prothrombin complex factors (II, VII, IX, X, Protein C and Protein S) are present at the level normally expected in plasma from patients under oral anticoagulant treatment. The PIVKA inhibitors (Proteins induced by Vitamin K Absence) are present as in patient's plasma. No preservatives are added to the lyophilised plasmas.

mL	Plasma	Description
1 x 1	Level 1	Lyophilized anticoagulated plasma, INR about 2.0
1 x 1	Level 2	Lyophilized anticoagulated plasma, INR about 3.0
1 x 1	Level 3	Lyophilized anticoagulated plasma, INR about 4.0

Materials required (not supplied with the kit)

- Pipettes
- Distilled water
- Calibrator and reagents
- REF 5010004 AK-Calibrant 4 x 1 mL
- REF 5003220 TECHNOCLOT® PT Plus 5 x 4 mL
- REF 5003222 TECHNOCLOT® PT Plus 10 x 10 mL
- REF 5005032 TECHNOCLOT® PT Owren manual 10 x 4 mL
- REF 5005037 TECHNOCLOT® PT Owren manual 10 x 10 mL
- REF 5005044 TECHNOCLOT® PT Owren automated 10 x 4 mL
- REF 5005046 TECHNOCLOT® PT Owren automated 10 x 10 mL

Warning and precautions

- IVD - for in vitro diagnostic use.
- This verification set is intended for use by personnel trained in laboratory procedures and universal precautions for use of chemicals and potentially biohazardous substances.
- All human blood or plasma products as well as samples must be considered as potentially infectious. They have to be handled with appropriate care and in strict observance of safety regulations. The rules pertaining to disposal are the same as applied to disposing hospital waste.
- Verification plasma pools are made from human blood and any individual plasma involved in the procedure is HbsAg, HIV 1/2 Ab and HCV-Ab-negative as tested by FDA approved methods. However, all human blood products should be handled as potentially infectious material.
- A Material Safety Data Sheet for this product is available from www.technoclone.com

Stability and storage

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8°C. Stability after reconstitution:

18...25°C	4°C	≤-20°C
4 hours	8 hours	1 month

The vials can be frozen only once.

TEST PROCEDURE

- All reagents have to reach room temperature (18...25°C) before mixing.
- Open the vial and reconstitute the contents in 1mL of distilled by carefully rotating the vial until the product is completely reconstituted (avoid frothing). Allow them to stand for 30 minutes at room temperature before use. Invert to mix before use (avoid foaming).

Reconstituted plasma should be treated in the same way as patient samples, using the instructions given for the test. Following the ISTH^{2,3} and CLSI⁴ guidelines run the AK-Controls in duplicate over a minimum of three days/run with you PT/instrument system.

LIMITATION OF THE TEST

A new verification is required for each batch of thromboplastin and for each instrument used. Also following a major service of either instrument or equipment, the system should be revalidated. The figures in the batch table are only applicable to the indicated lot number of AK-Verification Kit and the method has to adhere strictly as prescribed for the respective reagent. The AK-Verification Kit is only as good as the assay system. Possible problems with parts of the test system can be indicated by not achieving the reference values.

INTERPRETATION

Technoclone provides a free evaluation sheet on the website www.technoclone.com Following information need to be entered into the evaluation sheet:

- Thromboplastin and instrument information
- Lot specific or locally established ISI value (International Sensitivity Index),
- MNPT value (Mean Normal Prothrombin Time),
- INR reference values from the lot specific batch table.
- INR results of verification plasmas,

The evaluation sheet will automatically calculate the mean INR and CV for each level as well as the accuracy.

The CVs of duplicates and between runs need to be within 10% and accuracy within ±15% of reference values. The within-run precision was assessed over multiple runs using a specific lot of reagent and AK-Control plasmas. The coefficients of variation obtained in this study were <6% for PT [sec] for all levels. ISTH and CLSI documents state that the difference between duplicate PT [sec] as well as the difference between the mean of the run (day-to-day) results must be within 10%.

Standardisation

The INR's for AK-Verification Kit are determined according to WHO (WHO_TRS_889_A3_1999) guidelines by Technoclone and/or by external quality assessment.

LITERATURE

1. WHO Technical Report Series, No. 889, 1999, Annex 3, Guidelines for thromboplastins and plasma used to control oral anticoagulant therapy
 2. British Journal of Haematology, 1999, 104, 455-460, An international collaborative study on the INR calibration of freeze-dried reference plasmas. ANTHONY R. HUBBARD, SARAH M. L. MARGETTS, LYNNE J. WELLER, JO MACNAB* AND TREVOR W. BARRROWCLIFFE, Division of Haematology and *Informatics Laboratory, National Institute for Biological Standards and Control, Potters Bar
 3. J Thromb Haemost 2004; 2: 1946-53. Guidelines on preparation, certification, and use of certified plasmas for ISI calibration and INR determination. A. M. H. P. VAN DEN BESSELAAR, T. W. BARRROWCLIFFE*, L. L. HOUBOUYAN-RE* VEILLARD, J. JESPERSEN, M. JOHNSTON, S. L. POLLER and A. TRIPODI* ON BEHALF OF THE SUBCOMMITTEE ON CONTROL OF ANTICOAGULATION OF THE SCIENTIFIC AND STANDARDIZATION COMMITTEE OF THE ISTH
 4. H54-A Vol. 25 No. 23 2005, Procedures for Validation of INR and Local Calibration of PT/INR Systems: Approved Guideline
- For more information please contact Technoclone or your local distributor.

SYMBOL KEYS

	Manufacturer		CE-mark
	Storage temperature		In vitro diagnostic use
	Expiry date		Lot
	Consult instructions for use		Catalogue number
	Biological risk		Global Trade Item Number

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