Coagulation Control N

GB

IT

PT

SE

NO

BG

RU

HR

DE

ES

FR

DK




GR

CZ

SRB

REF 5020040 Coagulation Control N 5 x 1 mL

REF 5020050 Coagulation Control N 50 x 1 mL

Symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Klíčová slova / Značenje simbola			
	Manufacturer / Hersteller / fabbricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produsent / Κατασκευαστής / Производитель / Производител / Proizvođač		Expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace / срок годности / datum expirace / Rok trajanja
	Storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / opbevaringstemperatur / Oppbevaringstemperatur / Θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		Consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg bruksveiledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potfeba fidiť se instrukciami / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj uputstvo pre upotrebu
	CE-mark / CE-Kennzeichnung / marchio CE / marca de CE / Simbolo da CE / marquage CE / CE-märkning / CE-mærket / CE-merke / CE-σημάδι / CE марка / CE-označení / маркировка CE / značka CE / CE-marka		Determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelser / Bestemmelser / προσδιορισμοί / брой тестове / stanovení / определний / роčet stanovení / Definicija
	Distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / дестилирана вода / destilovaná voda / дистилюванна вода / destilovaná voda / Destilisaná Voda		Lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / партиа / партида номер / šarže / лот / šarže / Serija
	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionsbuffert / Reaktionsbuffer / Reaksjonsbuffer / διάλυμα αντίδρασης / Реакционный буфер / Рабочий буферный раствор / Reakční pufr / Reakcioni pufer		Microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροπιλοδότησης / Микротитерная плака / Микропланшет / Mikrotitrační destička / Mikrotitracione ploče
	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator		Catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat		Ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk / έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola		Stop solution / Stopplösning / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stoppløsning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп раствор / Стоп-раствор / Zastavovací roztok / Stop solucija
	Dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / späd eller upplös i / fortyndes eller opløses i / Fortyndes eller opløses i / αραιωση ή διάλυση σε / разворете или разрежете с / zředit anebo rozpustit v / разбавить или растворить в / nafedte nebo rozpustte v / razrediti ili rastvoriti u		Substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
	Incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationsbuffert / Inkubationsbuffer / Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационный буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer		Washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningsskoncentrat / Vaskeopløsningsskoncentrat / vaskeløsningsskoncentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промывочного раствора / Koncentrát promývacího roztoku / Koncetrat solucije za ispiranje
	In vitro diagnostic use / in vitro Diagnostikum / diagnostico in vitro / diagnóstico in vitro / diagnóstico in vitro / diagnostic in vitro / för in vitro diagnostik / in vitro diagnostik / in vitro diagnostisk bruk / ρηση διαγνωστικής εντόβ σωλήνα / за ин витро диагностика / pro in vitro diagnostiku / использовать для диагностики in vitro / diagnostický prostředek in vitro / in vitro diagnostika		

PRODUCT DESCRIPTION**INTENDEND USE**

For use as accuracy control of screening tests (PT/aPTT/TT/Fibrinogen) and determination of various coagulation factors within normal coagulation activity.

Precision control, i.e. the recording of the day to day variation, is an important procedure in coagulation testing; checking the accuracy of the test system is particularly important.

Coagulation Control N can be used for:

- Accuracy control of all screening tests, i.e. prothrombin time (PT), partial thromboplastin time (aPTT), thrombin time (TT), fibrinogen determination (Clauss), and of various coagulation factors in the normal coagulation activity.
- Precision control of all parameters listed in the table.

COMPOSITION

Coagulation Control N is prepared from selected citrated plasma donations of healthy donors. The clotting activity is normally distributed. Thereby the "average" presence of all coagulation factors and inhibitors is secured. Coagulation Control N contains stabilizers but no bactericide additives.

MATERIAL REQUIRED (not supplied with the kit)

- Pipettes
- Distilled water

WARNING AND PRECAUTIONS

- IVD for *in vitro* diagnostic use
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- Each single donor plasma and each lot of Coagulation Control N are tested and found negative for HbsAg, HIV 1/2 Ab and HCV Ab. However, universal precautions (treating all human source materials as if potentially infectious) should be exercised.

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C.

Stability after reconstitution:

RT* 4 hours	+2...8 °C 8 hours	-20 °C 1 month
----------------	----------------------	-------------------

The vials can be only frozen once.

Upon storage, caps should be screwed tightly.

*=room temperature ¹

TEST PROCEDURE**PREPARATION AND PERFORMANCE OF THE TEST**

- Open the vial carefully and reconstitute the contents in 1 mL of distilled water by carefully rotating the vial until the product is completely reconstituted (avoid frothing).
- Allow the reconstituted Coagulation Control N to stand for 10 minutes at room temperature before use ². Invert to mix before use (avoid foaming).
- Treat the reconstituted plasma as a citrated patient sample according to the instructions of the respective test.

ANALYSES RESULTS**EVALUATION OF THE RESULTS**

The figures in the table are only applicable to the indicated lot number of Coagulation Control N and the method has to adhere strictly as prescribed for the respective reagent.

STANDARDIZATION AND TRACEABILITY

The reported values were determined over multiple runs using a specific lot of reagent and against an Internal Reference Standard which is traceable to the current International Standards, identified in the acceptance range table.

For the tests where International Standards are not available, these parameters have been assigned against an Internal Reference Standard which is traceable to a frozen normal plasma pool of 100 donors.

LIMITATION OF THE TEST

The figures in the table are only applicable to the indicated lot number of Coagulation Control N and the method has to adhere strictly as prescribed for the respective reagent.

The guidelines for evaluation given in this leaflet have been established assuming that a constant sensitivity of different batches is guaranteed by the manufactures of each reagent.

LITERATURE

Please contact Technoclone or your local distributor.

¹ Use the reconstituted plasma in single factor determinations within 2 hours stored at room temperature.

² For standardization a reconstitution time of 30 minutes is recommended.

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