

**Your instrument of choice for classical haemostasis routine and speciality testing.  
No matter if your focus is on routine diagnostics or on research projects.**

Whether running patient samples in a routine setting or in deeply specialized research projects the new Ceveron c100 is the perfect addition to your lab. We at Technoclone managed to improve the ease of use of this a fast, robust and precise instrument with our new, touch screen optimised software.

In combination with our high quality and long on board stable Technoclone reagents, you can now achieve even more performance in a minimal timeframe.

With no more space needed than our first generation Ceveron alpha, we have refined our Ceveron 100 series platform into the perfect solution for your haemostasis speciality testing, without loosing focus on routine assays.

chromogenic assays  
turbidimetric assays  
clotting assays



# CEVERON c100

**Testing classical clotting, chromogenic and turbidimetric parameters without limitation.**

## **The perfect solution for routine and research analytics**

- Easy operating with the new, touch screen optimized software
- High capacity result storage including used reagent lots and raw data
- Easy export of result and raw data files
- Open system with simple test application programming

## **Unique robustness and precision**

- Very low cv in combination with Technoclone reagents
- High quality LED optics detect even smallest reaction changes in the cuvette
- Temperature monitoring of each individual measurement

## **Secure and high traceable**

- Different password protected security levels available
- All activities of the Ceveron c100 are tracked and stored in the audit trail
- Remote service access opportunity

- Global Tests PT, PT Owren, aPTT, Fibrinogen, Thrombin time
- Anticoagulant Management Apixaban, Argatroban, Arixtra, Dabigatran, Edoxaban, Orgaran, Rivaroxaban and Heparins
- Coagulopathy Factor Assays (incl. HMWK and Prekallikrein Plasma) FVIII Inhibitor, FVIII:C
- Thrombophilia Protein C, Protein S, AT, APC-R, C1-INH, Lupus
- Thrombosis D-Dimer

Product	Description	REF	Package
Ceveron c100	Fully automated coagulation analyzer for <b>clotting</b> , <b>chromogenic</b> and <b>turbidimetric assays</b>	9822010	1 pc.

CEV007E.002 11/2021

### **Technoclone Herstellung von Diagnostika und Arzneimitteln GmbH**

Brunner Str. 67 1230 Vienna Austria  
Phone: +43 1 86373-0 Fax: +43 1 86373-44  
products@technoclone.com www.technoclone.com



*Your partner in Coagulation !*

Ceveron is a registered trademark of Technoclone Herstellung von Diagnostika und Arzneimitteln GmbH.



NEW !



technoclone  
TC

## CEVERON 100 series

First fully automated coagulation analyzer series for clotting, chromogenic, turbidimetric assays, thrombin generation and quenching technology (FXIII & ADAMTS13 Activity)




*An explosion of possibilities !*

# The evolution of fully automated testing continues

**The first fully automated instrument series where you can combine classic haemostasis diagnostics with Thrombin Generation, Factor XIII Act\* and ADAMTS13 Act\***

Simultaneous measurement of five technologies.

Up to five technologies clotting, chromogenic, turbidimetric, thrombin generation and quenching can be run at the same time.

	Clotting	Chromogenic	Turbidimetric	Thrombin Generation	FXIII Act*	ADAMTS13 Act*
	✓	✓	✓			
	✓	✓	✓	✓		
	✓	✓	✓	✓	✓	✓



ADAMTS13 Act  
FXIII Act  
Thrombin Generation  
chromogenic assays  
turbidimetric assays  
clotting assays

\*Quenching Module




## High precision and flexibility

Rerun and reflex testing in combination with very low cv in routine and speciality testing, guarantee secure results for perfect diagnostics.



Having a minimal patient and reagent dead volume, adult's and children's samples can be run without any differentiation.

Due to the workflow of Ceveron 100 series an exact incubation time on each measurement is guaranteed.



Fibrinogen  
Fibrinogen ↗  
Fibrinogen ↘


Easy access to all raw data, help comparing patient results, QC results and calibration curves.



## Complete QC monitoring

- Westguard Rules
- Laboratory specific confidence range
- 3s range
- Lifetime of validation adjustable

Using values made available by Technoclone, via batch tables, local established results or pre-period results of new QC lots.



With a throughput of 75 PT/hour, a TGA result below 20 minutes and an ADAMTS13 **Act** result within 25 minutes, Ceveron 100 series instruments are the perfect solution for fast and precise diagnostic of patients.

Ceveron 100 series is ready to use within 5 minutes after switching on the instrument.

No daily maintenance needed. Just 15 minutes once per week.




## High traceability

Sample and reagent handling is fully traceable. With no difference in barcoded, or non barcoded sample and reagent handling.

All used reagent lots, for each measurement can be sent to LIS.

Actual temperatures during the pipetting, incubation and reaction of each sample are shown, and stored on Ceveron 100 series systems.



## System and result security

Three security levels ensure, that just qualified users can run samples, perform calibrations or change settings.


On each security level, additional administrator rights can be set.

Due to the unlimited sample result storage, no QC or patient result can be lost due to limited storage capacity.

All results are stored on the Ceveron, until they are deleted manually.

Automatic or manual backup including raw data and system status.

The integrated logbook tracks all actions, created by the instrument, a user or the service technician.



## Host connection

All Ceveron 100 series instruments are equipped with Windows 10 and use ASTM 1394 bidirectional host connection.

Obviously the Ceveron instruments can also operate in offline mode, as easy as connected to a LIS.


















## New intuitive icon based software

No matter which Ceveron 100 series instrument, they are using the same intuitive icon based, touch screen optimized software.


Operating the Ceveron 100 series instruments is even more simplified.



With our colour coded Ceveron symbol, the actual system status can be spotted at the first glimpse.



With just one click you can place samples, mark them as  Then  the reagents or just  their actual status. After  your cuvettes, just  your sample testing. Meanwhile you can check your  results or  your last calibration curves. After receiving the test results you can compare your, ,  the results manually or automatically to your LIS, or  them. If you have new lots of QC material, you can use the  or enter the target values manually. When performing a new calibration, after  you are  to run your samples. If you have written your own tests, you can easily them  to a second or third  instrument.

All patient, control and calibration results will be stored on the Ceveron until you  and  them.

All steps are as easy as above, for Quenching and  parameters.

Occasionally you have to  your system solution. Once a week you have to perform some  which claims 15 minutes of your valuable time.

Although we really appreciate all our Ceveron users, you will have a  visit just once per year.

Needless to say, we can also keep in  using  connection, if you need further

assistance, working with your new  instruments.

# CEVERON 100 series SPECIFICATIONS

Test Menu	Clotting, chromogenic, turbidimetric assays TGA Module*, Quenching Module*
Optics	4 channels with 4 wavelengths each wavelengths: 405, 540, 630, 740 nm 4 channels TGA Module* 4 channels Quenching Module*
Samples on board	36 (24 primary tubes +12 for sample cups)
Reagents on board	32 14 (12-18 °C), 10 (RT) 8 buffers and diluents
Cuvettes on board	84
Tests per sample	36
Sample volume	5 - 400 µL
Reagent volume	5 - 400 µL
Throughput	75 PT/h
STAT	Yes, up to 8 samples
Quality control	Integrated, user defined QC limits, Westgard rules
Patient data storage	unlimited
Operating system	Windows 10
Dimensions (analyzer)	750 mm x 675 mm x 420 mm (w x d x h)
Weight	40 kg

Sample barcode reader	✓
Reagent barcode reader	✓
Sample and reagent level sensor	✓
Continuous operation	✓
Sample predilution	✓
Calibration curve predilution	✓
Reaction curve availability	✓
Security system	✓
Primary and secondary vials	✓
Bi-directional interface	✓
Host query function	✓
Events log system	✓
Rerun and reflex testing	✓
User defined tests	✓

\* depending on instrument

Product	Description	REF	Package
<b>Ceveron c100</b>	Fully automated coagulation analyzer for <b>clotting, chromogenic</b> and <b>turbidimetric assays</b>	9822010	1 pc.
<b>Ceveron t100</b>	Fully automated coagulation analyzer for <b>clotting, chromogenic, turbidimetric and fluorogenic</b> assays (TGA Module)	9822110	1 pc.
<b>Ceveron s100</b>	Fully automated coagulation analyzer for <b>clotting, chromogenic, turbidimetric fluorogenic</b> and <b>quenching</b> assays (TGA Module and Quenching Module)	9822210	1 pc.

CEV006E.003 11/2021

**Technoclone Herstellung von Diagnostika und Arzneimitteln GmbH**  
 Brunner Str. 67 1230 Vienna Austria  
 Phone: +43 1 86373-0 Fax: +43 1 86373-44  
 products@technoclone.com www.technoclone.com



*Your partner in Coagulation!*

Ceveron is a registered trademark of Technoclone Herstellung von Diagnostika und Arzneimitteln GmbH.