



TECHNOZYM® vWF:CBA ELISA

Calibrators and Controls

LOT RB45C00

3015987-178

Table • Tabelle • tabella • tabla • tabela • tabell

DIL 0.5mL AQUA

REF CALIBRATORS	LOT		Sollwert/ Nominal value	Kalibrierung durch®, nach Calibration by®, acc, to	
vWF CAL 1 TH10001ND	1ND32C0	2025-12-31	1.68 IU/mL	IS 07/316	
vWF CAL 2 TH10002ND	2ND31C0	2025-12-31	0.83 IU/mL	IS 07/316	
vWF CAL 3 TH10003ND	3ND31C0	2025-12-31	0.39 IU/mL	IS 07/316	
vWF CAL 4 TH10004ND	4ND31C0	2025-12-31	0.19 IU/mL	IS 07/316	
vWF CAL 5 TH10005ND	5ND31C0	2025-11-30	0.00 IU/mL	IS 07/316	

REF CONTROLS	LOT		Sollwert/ Nominal value	Vertrauensbereich Confidence range	Kalibrierung durch®, nach Calibration by®, acc, to
vWF CONT H TH1000AND	AND32C0	2025-11-30	0.96 IU/mL	0.72 – 1.20 IU/mL	IS 07/316
vWF CONT L TH1000BND	BND31C0	2025-11-30	0.41 IU/mL	0.25 – 0.57 IU/mL	IS 07/316

