




# Technoview Organan CAL

**LOT** 0131C00

 **2025-01-31**

3015175- 019 19/01/2023

Test test	Lot	Einheit unit	Richtwert (Vertrauensbereich) expected value (confidence range)	Std.
Technochrom Anti Xa	CAL 1 5L11C01	IU/mL	0.00	Internal Reference Standard
	CAL 2 5M12C01	IU/mL	0.36	
	CAL 3 5N12C01	IU/mL	0.70	
	CAL 4 5O12C01	IU/mL	1.24	
	CAL 5 5P12C01	IU/mL	1.84	