

Mit 0.5 mL Aqua dest. rekonstituieren / reconst. with 0.5 mL Aqua dest.

| Artikel / Article CALIBRATORS | Lot / batch | Ablaufdatum / expiry date | Sollwert / Nominal value Adamts - 13 Activity | Sollwert / Nominal value Adamts - 13 Antigen |
|--|-------------|---------------------------|--|---|
| Adamts - 13 Calibrator 1 Art.Nr.: 5450610 | 1RT81B1 | 2021-01-31 | 0.79 IU/mL | 0.88 IU/mL |
| Adamts - 13 Calibrator 2 Art.Nr.: 5450612 | 2RT81B1 | 2021-01-31 | 0.41 IU/mL | 0.48 IU/mL |
| Adamts - 13 Calibrator 3 Art.Nr.: 5450614 | 3RT81B1 | 2021-01-31 | 0.20 IU/mL | 0.23 IU/mL |
| Adamts - 13 Calibrator 4 Art.Nr.: 5450616 | 4RT81B1 | 2021-01-31 | 0.10 IU/mL | 0.12 IU/mL |
| Adamts - 13 Calibrator 5 Art.Nr.: 5450618 | 5RT81B1 | 2021-01-31 | 0.00 IU/mL | 0.00 IU/mL |

| Artikel / Article CONTROLS | Lot / batch | Ablaufdatum / expiry date | Sollwert / Nominal value | | Vertrauensbereich Confidence range |
|---|-------------|---------------------------|--------------------------|------------|---------------------------------------|
| Adamts - 13 High Control Art.Nr.: 5450620 | RRY81B1 | 2021-01-31 | Antigen | 0.71 IU/mL | 0.53 – 0.89 IU/mL |
| | | | Activity | 0.59 IU/mL | 0.44 – 0.74 IU/mL |
| Adamts - 13 Low Control Art.Nr.: 5450622 | RRZ81B1 | 2021-01-31 | Antigen | 0.18 IU/mL | 0.11 – 0.25 IU/mL |
| | | | Activity | 0.16 IU/mL | 0.10 – 0.22 IU/mL |

IU/mL wurden gegen den WHO 1st international standard ADAMTS-13 plasma (12/252) kalibriert.
IU/mL have been calibrated against the WHO 1st international standard ADAMTS-13 plasma (12/252)