

# Ceveron TGA RC High Kit

LOT **B543C00**

 **2025-08-31**

3013832- 021 04/10/2024

Ceveron 100 series barcode



Kontrolle control	Lot	Parameter	Einheit unit	Richtwert (Vertrauensbereich) expected value (confidence range)	Std.
Ceveron TGA Control High	1Y42C00	Peak thrombin	nM	469 ( 352 - 586 )	IS 19/188
Ceveron TGA Control Low	2742C00			272 ( 204 - 340 )	IS 19/188