



**Technoview Organan CAL**
**LOT 0121C00**
 **2023-08-31**
3015175- 018 12/01/2022

| Test test           | Lot           | Einheit unit | Richtwert (Vertrauensbereich) expected value (confidence range) | Std.                        |
|---------------------|---------------|--------------|---|-----------------------------|
| Technochrom Anti Xa | CAL 1 5L11C00 | IU/mL        | 0.00  | Internal Reference Standard |
|                     | CAL 2 5M12C00 | IU/mL        | 0.36  |                             |
|                     | CAL 3 5N12C00 | IU/mL        | 0.70  |                             |
|                     | CAL 4 5O12C00 | IU/mL        | 1.24  |                             |
|                     | CAL 5 5P12C00 | IU/mL        | 1.84  |                             |