

technoclone  
TC

# FACTOR DEFICIENT PLASMA


**Faktor II, V, VII, VIII, IX, X, XI and XII**

**Fitzgerald and Fletcher Trait Plasma**

*A clear view into the future!*

# FACTOR ASSAYS

Factors of the intrinsic Pathway are tested in one-stage clotting assays with activated Partial Thromboplastin Time (aPTT) reagents using plasma deficient of Factor VIII, IX, XI and XII as well as Fitzgerald and Fletcher Trait Plasma.



Factors of the extrinsic and Common Pathway are tested in one-stage clotting assays with Prothrombin Time (PT) reagents using plasma deficient of Factor II, V, VII and X.

Product	Description	REF	Package
<b>Factor II deficient plasma</b>	Lyophilized Factor II deficient plasma, immunodepleted	5114008	5 x 1 ml
<b>Factor V deficient plasma</b>	Lyophilized Factor V deficient plasma, immunodepleted	5134004	5 x 1 ml
<b>Factor VII deficient plasma</b>	Lyophilized Factor VII deficient plasma, immunodepleted	5144015	5 x 1 ml
<b>Factor X deficient plasma</b>	Lyophilized Factor X deficient plasma, immunodepleted	5174006	5 x 1 ml
	Lyophilized Factor X deficient plasma, native	5174004	5 x 1 ml
<b>Factor VIII deficient plasma</b>	Lyophilized Factor VIII deficient plasma, immunodepleted	5154002 5154004	5 x 1 ml 50 x 1 ml
	Lyophilized Factor VIII deficient plasma, native	5154007 5154016	5 x 1 ml 50 x 1 ml
<b>Factor IX deficient plasma</b>	Lyophilized Factor IX deficient plasma, immunodepleted	5164003 5164004	5 x 1 ml 50 x 1 ml
	Lyophilized Factor IX deficient plasma, native	5164008 5164016	5 x 1 ml 50 x 1 ml
<b>Factor XI deficient plasma</b>	Lyophilized Factor XI deficient plasma, native	5184004	5 x 1 ml
	Lyophilized Factor XI deficient plasma, immunodepleted	5184006	5 x 1 ml
<b>Factor XII deficient plasma</b>	Lyophilized Factor XII deficient plasma, native	5194008	5 x 1 ml
	Lyophilized Factor XII deficient plasma, immunodepleted	5194007	5 x 1 ml
<b>Fitzgerald Trait Plasma</b>	Lyophilized high molecular weight kininogen deficient plasma, immunodepleted	5204006	2 x 1 ml
<b>Fletcher Trait Plasma</b>	Lyophilized prekallikrein deficient plasma, native	5205006	2 x 1 ml
Auxiliary reagents required and not supplied with the kit			
<b>Citrate Sodium Chloride Buffer</b>	Dilution buffer for determination of coagulation factors II, V, VII and X with PT reagents	5400045 5400047	60 ml 25 ml
<b>Imidazole Buffer</b>	Dilution buffer for use in factors VIII, IX, XI, XII and factor VIII Inhibitor, Protein C and Fibrinogen.	5410008 5410007	25 ml 9 x 25 ml
<b>CaCl<sub>2</sub> 50 mM</b>	Dilution buffer for determination of Fitzgerald and Fletcher Trait Plasma in clotting method.	5279027	25 ml

