

Safety Data Sheet



Temp.: TF0366.01
Att 01 to DC0068

Technoview Edoxaban CON L

SDS No. P700.02
p 1/6 EN

ACCORDING TO EC-REGULATIONS 1907/2006 (REACH), 1272/2008 (CLP) & 2020/878

SECTION 1: IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

1.1 Product identifier

Product Name/Reference No. Technoview Edoxaban CON L 5 x 1 mL REF 5090251

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified Use(s) For professional users only
Uses Advised Against Not known.

1.3 Details of the supplier of the safety data sheet

Manufacturer

Company Identification Technoclone Herstellung von Diagnostika und Arzneimitteln GmbH
Address of Manufacturer Brunner Strasse 67
1230 Vienna
Austria
Postal code 1230
Telephone: +43 1 86373-0
Fax +43 1 86373-44
E-mail products@technoclone.com
Office hours 8:00 – 16:00

1.4 Emergency telephone number

Emergency Phone No. +43 1 86373-0
Contact products@technoclone.com
National response center
Address NHS Direct
Emergency Phone No. +44 111

SECTION 2: HAZARDS IDENTIFICATION

2.1 Classification of the substance or mixture

Regulation (EC) No. 1272/2008 (CLP) Not classified as dangerous for supply/use

2.2 Label elements

Product Name According to Regulation (EC) No. 1272/2008 (CLP)
Technoview Edoxaban CON L
Hazard Pictogram(s) None.

2.3 Other hazards

None known.

2.4 Additional Information

Contains human derived material - potentially biohazardous

SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

3.1 Substances

Not applicable.

3.2 Mixtures

Product does not contain any hazardous ingredients.

SECTION 4: FIRST AID MEASURES

4.1 Description of first aid measures

Inhalation If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing.
Skin Contact Wash skin with water.
Eye Contact Flush eyes with water for at least 15 minutes.

Safety Data Sheet



SDS No. P700.02
p 2/6 EN

Temp.: TF0366.01
Att 01 to DC0068

Technoview Edoxaban CON L

Ingestion Wash out mouth with water.

4.2 Most important symptoms and effects, both acute and delayed

None anticipated. Treat symptomatically.

4.3 Indication of any immediate medical attention and special treatment needed

Unlikely to be required but if necessary treat symptomatically.

SECTION 5: FIREFIGHTING MEASURES

5.1 Extinguishing media

Suitable Extinguishing media As appropriate for surrounding fire.
Unsuitable extinguishing media None.

5.2 Special hazards arising from the substance or mixture

None anticipated. Heating may cause decomposition.

5.3 Advice for firefighters

As appropriate for surrounding fire.

SECTION 6: ACCIDENTAL RELEASE MEASURES

6.1 Personal precautions, protective equipment and emergency procedures

Provide adequate ventilation. Wear suitable gloves if prolonged skin contact is likely.

6.2 Environmental precautions

Do not release large quantities into the surface water or into drains.

6.3 Methods and material for containment and cleaning up

Sweep spilled substances into containers if appropriate moisten first to prevent dusting.

6.4 Reference to other sections

See also Section 8, 13.

SECTION 7: HANDLING AND STORAGE

7.1 Precautions for safe handling

Not known.

7.2 Conditions for safe storage, including any incompatibilities

Storage temperature +2°C to +8°C
Storage life Stable under normal conditions.
Incompatible materials None known.

7.3 Specific end use(s)

In vitro diagnostic use

SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

8.1 Control parameters

8.1.1 Occupational Exposure Limits No Occupational Exposure Limit assigned.

8.2 Exposure controls

8.2.1. Appropriate engineering controls Ensure adequate ventilation.

8.2.2. Personal protection equipment

Eye Protection Wear eye protection with side protection (EN166).
Skin protection Gloves recommended
Respiratory protection Normally no personal respiratory protection is necessary.
Thermal hazards None known.

Safety Data Sheet



Temp.: TF0366.01
Att 01 to DC0068

Technoview Edoxaban CON L

SDS No. P700.02

p 3/6 EN

8.2.3. Environmental Exposure Controls Do not release large quantities into the surface water or into drains.

SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

9.1 Information on basic physical and chemical properties

| | |
|----------------------------------------------------------|--------------------------------------------------------------------------|
| Physical state | solid (lyoph.) |
| Colour | white-yellowish |
| Odour | odourless |
| Melting point/freezing point | Not known. |
| Boiling point or initial boiling point and boiling range | Not applicable. |
| Flammability | Not known. |
| Lower and upper explosion limit | Not known. |
| Flash Point | Not applicable. |
| Auto-ignition temperature | Not known. |
| Decomposition Temperature | Not known. |
| pH | Not known. |
| Kinematic Viscosity | Not known. |
| Solubility | Solubility (Water) : readily miscible. Solubility (Other) : Not known |
| Partition coefficient n-octanol/water (log value) | Not known. |
| Vapour pressure | Not known. |
| Density and/or relative density | Not known. |
| Relative vapour density | Not known. |
| Particle characteristics | Not known. |

9.2 Other information

None.

SECTION 10: STABILITY AND REACTIVITY

10.1 Reactivity

None anticipated.

10.2 Chemical Stability

Stable under normal conditions.

10.3 Possibility of hazardous reactions

No hazardous reactions known if used for its intended purpose.

10.4 Conditions to avoid

None anticipated.

10.5 Incompatible materials

Not known.

10.6 Hazardous decomposition products

No hazardous decomposition products known.

SECTION 11: TOXICOLOGICAL INFORMATION

11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008

| | |
|--------------------------------|--------------------|
| Acute toxicity - Ingestion | Low acute toxicity |
| Acute toxicity - Skin Contact | Low acute toxicity |
| Acute toxicity - Inhalation | Low acute toxicity |
| Skin corrosion/irritation | Non-irritant |
| Serious eye damage/irritation | Non-irritant |
| Skin sensitization data | Not classified |
| Respiratory sensitization data | Not classified |

Safety Data Sheet



Temp.: TFF0366_01
Att 01 to D00068

Technoview Edoxaban CON L

SDS No. P700.02
p 4/6 EN

| | |
|--------------------------|----------------------------------------------|
| Germ cell mutagenicity | There is no evidence of mutagenic potential. |
| Carcinogenicity | No evidence of carcinogenicity. |
| Reproductive toxicity | Not classified |
| Lactation | Not classified |
| STOT - single exposure | Not classified |
| STOT - repeated exposure | Not classified |
| Aspiration hazard | Not classified |

11.2 Information on other hazards

Not known.

SECTION 12: ECOLOGICAL INFORMATION

12.1 Toxicity

| | |
|------------------------------------|--------------------------------|
| Toxicity - Aquatic invertebrates | Low toxicity to invertebrates. |
| Toxicity - Fish | Low toxicity to fish. |
| Toxicity - Algae | Low toxicity to algae. |
| Toxicity - Sediment Compartment | Not classified. |
| Toxicity - Terrestrial Compartment | Not classified. |

12.2 Persistence and degradability

Not known.

12.3 Bioaccumulative potential

Not known.

12.4 Mobility in soil

Not known.

12.5 Results of PBT and vPvB assessment

Not known.

12.6 Endocrine disrupting properties

None known.

12.7 Other adverse effects

Not known.

SECTION 13: DISPOSAL CONSIDERATIONS

13.1 Waste treatment methods

Dispose at suitable refuse site.

13.2 Additional Information

No special precautions are required for this product.

SECTION 14: TRANSPORT INFORMATION

Not classified as hazardous for transport.

14.1 UN number or ID number

Not applicable

14.2 UN proper shipping name

Not applicable

14.3 Transport hazard class(es)

Not applicable

14.4 Packing group

Not applicable

Safety Data Sheet

Technoview Edoxaban CON L



SDS No. P700.02
p 5/6 EN

Temp.: TF0366.01
Att 01 to DC0068

14.5 Environmental hazards

Not classified as a Marine Pollutant.

14.6 Special precautions for user

Not known

14.7 Maritime transport in bulk according to IMO instruments

Not known

SECTION 15: REGULATORY INFORMATION

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

European Regulations - Authorisations and/or Restrictions On Use

Candidate List of Substances of Very High Concern for Authorisation Not listed

REACH: ANNEX XIV list of substances subject to authorisation Not listed

REACH: Annex XVII Restrictions on the manufacture, placing on the market and use of certain dangerous substances, mixtures and articles Not listed

Community Rolling Action Plan (CoRAP) Not listed

Regulation (EC) N° 850/2004 of the European Parliament and of the Council on persistent organic pollutants Not listed

Regulation (EC) N° 1005/2009 on substances that deplete the ozone layer Not listed

Regulation (EU) N° 649/2012 of the European Parliament and of the Council concerning the export and import of hazardous chemicals Not listed

National regulations

Other Not known.

15.2 Chemical Safety Assessment

A REACH chemical safety assessment has not been carried out.

SECTION 16: OTHER INFORMATION

Date of Issue 2023-09-22 (replaces 358 v01 from 2017-05-10)
Date of Revision 2024-07-25, 16. reference to previous SDS no. added

LEGEND

Hazard Pictogram(s) None.
Precautionary Statement(s) None.

Acronyms
ATE : Acute Toxicity Estimate
CAS : Chemical Abstracts Service
CLP : Regulation (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures
DNEL : Derived No Effect Level
EC : European Community
EINECS : European Inventory of Existing Commercial Chemical Substances
LTEL : Long term exposure limit
PBT : Persistent, Bioaccumulative and Toxic
PNEC : Predicted No Effect Concentration

Safety Data Sheet



Temp.: TFF0366.01
Att 01 to D00068

Technoview Edoxaban CON L

SDS No. P700.02

p 6/6 EN

REACH : Registration, Evaluation, Authorisation and Restriction of Chemicals

STEL : Short term exposure limit

STOT : Specific Target Organ Toxicity

vPvB : very Persistent and very Bioaccumulative

Key literature references and sources Regulation (EC) No. 1272/2008 (CLP)
for data used to compile the SDS

Disclaimers

Information contained in this publication or as otherwise supplied to users is believed to be accurate and is given in good faith, but it is for the Users to satisfy themselves of the suitability of the product for their own particular purpose. Technoclone Herstellung von Diagnostika und Arzneimitteln GmbH gives no warranty as to the fitness of the product for any particular purpose and any implied warranty or condition (statutory or otherwise) is excluded except to the extent that exclusion is prevented by law. Technoclone Herstellung von Diagnostika und Arzneimitteln GmbH accepts no liability for loss or damage (other than that arising from death or personal injury caused by defective product, if proved), resulting from reliance on this information. Freedom under Patents, Copyright and Designs cannot be assumed.