

Mit 0.5 mL Aqua dest. rekonstituieren / reconst. with 0.5 mL Aqua dest.

Artikel/Article CALIBRATORS	Lot/batch	Ablaufdatum /expiry date	Sollwert / Nominal value Adamts-13 Antigen
Adamts - 13 Calibrator 1 Art.Nr.: 5450610	1RT61B0	2018-11-30	0.83 IU/mL
Adamts - 13 Calibrator 2 Art.Nr.: 5450612	2RT61B0	2018-11-30	0.48 IU/mL
Adamts - 13 Calibrator 3 Art.Nr.: 5450614	3RT61B0	2018-11-30	0.23 IU/mL
Adamts - 13 Calibrator 4 Art.Nr.: 5450616	4RT61B0	2018-11-30	0.13 IU/mL
Adamts - 13 Calibrator 5 Art.Nr.: 5450618	5RT61B0	2018-12-31	0.00 IU/mL

Artikel/Article CONTROLS	Lot/batch	Ablaufdatum/expiry date	Sollwert / Nominal value	Vertrauensbereich Confidence range
Adamts-13 Antigen High Control Art.Nr.: 5450620	RRY61B0	2018-12-31	0.63 IU/mL	0.47 – 0.79 IU/mL
Adamts-13 Antigen Low Control Art.Nr.: 5450622	RRZ61B0	2018-12-31	0.21 IU/mL	0.13 – 0.29 IU/mL

IU/mL wurden gegen den WHO 1st international standard ADAMTS-13 plasma (12/252) kalibriert.
 IU/mL have been calibrated against the WHO 1st international standard ADAMTS-13 plasma (12/252)

Qualitätskontrolle
 quality control




Technoclone Herstellung von
 Diagnostika und Arzneimitteln GmbH
 Brunner Str. 67 - 1230 Vienna, Austria